
MW Lens

CN202222833773.0
CN202223038253.7


Product Advantages

- Micro eccentric ring technique: Each ring is in a non concentric state to continuously provide differentiated stimulation to the periphery of the retina.
- Small central optical zone: Under the premise of ensuring clear vision in the central macular area, defocusing stimulation is placed as close as possible to the sensitive area of the retina to enhance the effect of myopia intervention.
- Optimization of ring spacing: The ring spacing in sensitive areas of the retina is relatively small, further enhancing the effect of myopia control.

Most Suitable Wearers

- Children and adolescents with progressive myopia
- Children and adolescents with insufficient reserve of farsightedness

Fundamental Hypothesis

By projecting the myopic defocus signal onto the periphery of the retina through a micro lens with positive additional refractive power, the axial growth of the eye is suppressed.

Micro eccentric ring technology

The use of micro eccentric ring technology ensures that the centers of each ring are not at the same point, providing differentiated stimulation for the retina periphery and avoiding the phenomenon of "static defocus saturation" caused by constant defocus within the same field of view, in order to provide sustained and long-term myopia intervention effectiveness.

Patents Provide Technical Protection

CN202222833773.0

CN202223038253.7

Clinical Trials Provide Scientific Support

Cooperative Medical Care:

Sinopharm Dongfeng General Hospital.

15-month clinical data results:

Compared with low scattering single vision lens, MW lens have a delayed refractive power increment of 0.53D, a delayed axial increment of 0.24mm, an effective control rate of 50% for equivalent spherical power, and an effective control rate of 56% for axial length.

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Stimulus intensity evaluation---DPME

DPME is used to evaluate the performance of lenslet lens based on defocus amount, filling rate, diameter, surface shape, and optical zone diameter of the lens.

product	DIMS	HAL	CAREA	CAREB	MW
DPME	8.124	11.140	4.553	7.604	12.562

Production Information

refractive index	1.59 (PC)	Lens diameter: • SE less than -8.00D, 70mm • SE is -8.25D~-9.50D, 65mm • SE is -9.75D~11.00D, 60mm
DS	-11.00D~+2.00D	
DC	-4.00D~0.00D	
SE	-11.00D	
coating	green film, super water film, anti blue light film	